Date: 15.02.2021 Bearbeiter: Harald Rösner Work area: Transmission

Electron Microscopy Workplace: Room E.05/07

## Safety Instruction **Titan Themis G3**

WWU Münster Nanofabrication Facility (MNF) CeNTech II Heisenbergstr. 11 48149 Münster

	Range of Application	
	Microscope Operation	
	Hazard for People and the environment	
	<ul> <li>High voltage: Danger due to high tension.</li> <li>Process gases: Mortal danger when inhaled. Suffocating in high concentrations.</li> <li>Gas bottles: Explosion hazard.</li> <li>Iionizing radiation: Full protection device. Info: Dr. H. Rösner</li> <li>Solvent: Chemical burn hazard.</li> <li>Cold liquid gas: Liquid nitrogen. Chokedamp. Burn hazard.</li> <li>Risk of falling: ramp</li> </ul>	Danger Liquid nitrogen
	Safety Measures and Rules of Behavior	
	<ul> <li>Operating the microscope is only allowed after previous instructions to the Titan Themis by Dr. Ostendorp and/or Dr. Rösner.</li> <li>Wear protective equipment (safety googles and gloves) when using liquid nitrogen.</li> </ul>	
	In case of any unclear issues contact the responsible person Dr. H. Rösner.	
	Action in Failure	
	In case of any malfunction or abnormality contact Dr. Ostendorp and Dr Rösner.	
	Action in Accident, First Aid	
-	<ul> <li>After contact with the eye: wash out for 10 minutes with floating water and open lid slit, immediately consult an eye doctor. After inhalation: get fresh air, immediately consult a doctor.</li> <li>Liquid nitrogen accident: Take off soaked clothes. In case of heavy burns contact emergency doctor.</li> <li>Inform appointed first aider.</li> <li>Immediately treat small injuries.</li> <li>Enter data in the accident book.</li> <li>In case of severe injuries call an emergency doctor (phone 112).</li> <li>Inform your superior.</li> </ul>	
	Maintenance and Disposal	
	Maintenance work may only be carried out by qualified and instructed staff.	